

**PART# 23046**



**Manufactured by:**



COMPASSIONATE SOLUTIONS

# Your First Few Days

As a result of using pads or condom catheters, your skin currently has **A LOT** of excess moisture. Fortunately for you, Men's Liberty™ just **LOVES** moisture, so it will soak it all up, until your skin returns to a normal moisture balance.

It can take up to three consecutive applications of Men's Liberty for you to get a full 24-hour wear time. **The most important thing** is to reapply another Men's Liberty right away – don't go back to other products as they will reintroduce all that moisture!

 **Call our nurses at (888) 920-9363 with any questions**

## Before Applying Men's Liberty

-  Begin by washing your anatomy with Ivory or Gold Bar Dial soap.
-  Make sure you are starting with skin that is completely clean and dry.
-  If you are uncircumcised or leak, use a KindKlamp™ penile clamp from BioDerm® to aid in application.
-  Seat yourself somewhere comfortable, with your Men's Liberty within easy reach.

**1** Begin by washing the head of the penis with Ivory or Gold Bar Dial soap. Allow skin to dry completely.



**2** Wipe the head of the penis with the BioPlus+™ wipe. Allow 15-20 seconds to dry until the skin is no longer tacky.

**3** Remove release papers and attach the faceplate to the tip of the penis with the center opening directly in line with the urinary opening. Smooth down petals.



**4** Apply the faceplate seal around the head of the penis flat and wrinkle free. Be sure the hydrocolloid is pressed into all the ridges and grooves in the skin. Press and hold for 10-15 seconds to secure the seal.



**For overnight use:** Unplug the bottom of the Men's Liberty pouch and connect your bedside bag. Hang the bed bag from your bed frame or rest against the floor. Keep below the bladder to ensure drainage. Do not twist or kink the tube.

*If irritation occurs, discontinue use and consult your healthcare professional.*

## Removing Men's Liberty

FreeDerm™ is the recommended remover for Men's Liberty. If you do not wish to use FreeDerm, Men's Liberty can be removed by soaking the seal in warm water or a wet washcloth until it is milky white and rolls gently from the skin.

*Never pull or tear from the skin.*



**BioDerm, Inc.**  
12320 73rd Court North  
Largo, FL, 33773, USA  
Orders: 800-373-7006  
Support: 800-864-6755  
CustomerCare@BioDermInc.com



**Authorized Representative:**  
Advena Ltd.  
Pure Offices, Plato Close,  
Warwick, CV34 6WE, UK

©2017 BioDerm, Inc. All rights reserved. BioDerm® is a registered trademark of BioDerm, Inc. Men's Liberty™, FreeDerm™, BioPlus+™ and KindKlamp™ are trademarks of BioDerm, Inc. Men's Liberty US patent #8,551,062 LT 90145 0

**NOT MADE WITH  
NATURAL RUBBER LATEX**